LIMITED WARRANTY

The Fingertip Pulse Oximeter with Wireless Connectivity has a 2 year warranty to be free of manufacturing defects in materials and workmanship under normal applications of the original owner. If this product becomes inoperable due to defect, return the product with all component pieces and proof of purchase to the address listed below.
This warranty does not cover any shipping/transport costs. This warranty does not apply if the product is subject to misuse, neglect, rough handling or damage.

Ship the unit prepaid and insured (at owner's option) to:
Thermor Ltd.
Attn: Returns Department
16975 Leslie Street
Newmarket, ON L3Y 9A1
www.biosmedical.com
Email: support@biosmedical.com

Limitations
This warranty covers all defects encountered in normal use of the equipment and does not apply in the following cases.
1) If the equipment has been serviced by anywhere other than a certified center.
2) Damage to the equipment due to mishandling, abuse, accident, or not following operating instructions.
3) INSTRUMENTS: Warranty does not extend to display face, batteries, carrying case, or lanyard.

For troubleshooting assistance, please contact Customer Service.
Our representatives are trained to provide assistance over the phone and may be able to resolve your problem without the need for a return. For the best possible support please have your unit readily available when calling us.
Service Phone: 1 (866) 517-2970
www.biosmedical.com

Appendix I

State	Prompt condition delay	Prompt signal generation delay
Low voltage prompt	1s	20ms
SpO_2 prompt	330ms	20ms
Pulse rate prompt	330ms	20ms
Probe error prompt	16ms	20ms

Appendix II EMC Guidance and Manufacturer Declaration

Table 1:

Guidance and manufacturer's declaration –electromagnetic emission	
The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The purchaser or the user of the device should assure that it is used in such environment.	
Emission test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B

Table 2:

Guidance and manufacturer's declaration-electromagnetic immunity		
The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The purchaser or the user of the Pulse Oximeter should assure that it is used in such environment.		
Immunity test	IEC60601 test level	Compliance level
Electrostatic discharge (ESD)	±8kV contact	±8kV contact
IEC 61000-4-2	± 15 kV air	±15kV air
Power frequency (50 / 60Hz) magnetic field IEC 61000-4-8	30 A/m	30A/m

Table 3:

Guidance and manufacturer's declaration – electromagnetic immunity		
The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer the user of the Pulse Oximeter should assure that it is used in such environment.		
Immunity test	IEC 60601 test level	Compliance level
Radiated RF IEC61000-4-3	10 V/m 80 MHz- 2.7 GHz	10 V/m80 MHz- 2.7 GHz

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Pulse Oximeter is used exceeds the applicable RF compliance level above, the Pulse Oximeter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Pulse Oximeter.
b) Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3V/m.

Table 4:

Guidance and manufacturer's declaration - electromagnetic Immunity							
The [Code SI] is intended for use in the electromagnetic environment specified below. The customer or the user of the Pulse Oximeter should assure that it is used in such an environment							
	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380 –390	TETRA 400	Pulse modulation b) 18 Hz	1,8	0,3	27
	450	430 –470	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0,3	28
	710	704 –787	LTE Band 13,17	Pulse modulation b) 217 Hz	0,2	0,3	9
	745						
	780						
	810	800 –960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0,3	28
	870						
	930						
	1720	1700 –1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3,4, 25; UMTS	Pulse modulation b) 217 Hz	2	0,3	28
	1845						
	1970						
	2450	2400 –2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0,3	28
	5240	5100 –5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0,2	0,3	9
	5500						
	5785						

NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.
b) The carrier shall be modulated using a 50 % duty cycle square wave signal.
c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation: $E = \frac{6}{d}\sqrt{P}$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

OPERATOR'S MANUAL

113PO

General Description

Oxygen binds to hemoglobin in red blood cells when moving through the lungs. It is transported throughout the body as arterial blood. A pulse oximeter uses two frequencies of light (red and infrared) to determine the percentage (%) of hemoglobin in the blood that is saturated with oxygen. The percentage is called blood oxygen saturation, or SpO_2. A pulse oximeter also measures and displays the pulse rate at the same time it measures the SpO_2 level.
Diagram of Operation Principle
1. Red and Infrared-ray Emission Tube
2. Red and Infrared-ray Receipt Tube

Precautions for Use

1. Before use, carefully read the manual.
2. Operation of the fingertip pulse oximeter may be affected by the use of an electrosurgical unit (ESU).
3. The fingertip pulse oximeter must be able to measure the pulse properly to obtain an accurate SpO_2 measurement. Verify that nothing is hindering the pulse measurement before relying on the SpO_2 measurement.
4. Explosive Hazard – **DO NOT** use this device in the environment of inflammable gases such as anesthetic.
5. **DO NOT** use the fingertip pulse oximeter in an MRI or CT environment.
6. This device is not intended for clinical diagnosis. If you are concerned with any results, review your measurements and associated symptoms with a licensed doctor or healthcare professional.
7. **DO NOT** use the fingertip pulse oximeter in situations where alarms are required. The device has no alarms. It is not for continuous monitoring.
8. **DO NOT** use the fingertip pulse oximeter in an explosive atmosphere.
9. The fingertip pulse oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.
10. In order to ensure correct sensor alignment and skin integrity, the maximum application time at a single site for our device should be less than half an hour.
11. When the device is on, **DO NOT** stare at the red and infrared light emitter (infrared light is invisible). This may be harmful to the eyes.
12. This device contains silicone, PCS, TPU, TPE, and ABS materials, whose biocompatibility has been tested in accordance with the requirements in ISO 10993-1, and has passed. Those allergic to silicone, PVC, TPU, TPE, or ABS should not use this device.
13. **DO NOT** sterilize the device using autoclaving, ethylene oxide sterilizing, or immersing the device in liquid. The device is not intended for sterilization.
14. Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
15. Portable and mobile RF communications equipment can affect medical electrical equipment.
16. This equipment is not intended for use during patient transport outside the healthcare facility.
17. This equipment should not be used adjacent to or stacked with other equipment.
18. It may be unsafe to:
 —use accessories, detachable parts and materials not described in the instructions for use
 —interconnect this equipment with other equipment not described in the instructions for use
 —disassemble, repair or modify the equipment
19. When the signal is not stable, the reading may be inaccurate. Please **DO NOT** reference.

Inaccurate Measurements May Be Caused By

1. Significant levels of dysfunctional hemoglobin (such as carbonyl - hemoglobin or methemoglobin).
2. Intravascular dyes such as indocyanine green or methylene blue.
3. High ambient light. Shield the sensor area if necessary.
4. Excessive patient movement.
5. High-frequency electrosurgical interference and defibrillators.
6. Venous pulsations.
7. Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.
8. The patient has hypotension, severe vasoconstriction, severe anemia, or hypothermia.
9. The patient is in cardiac arrest or is in shock.
10. Fingernail polish or false fingernails.
11. Weak pulse quality (low perfusion).
12. Low hemoglobin.

Product Features

1. SpO_2, pulse rate, and waveform values displayed
2. Automatic display direction
3. Free MediLink app to track measurements
4. Standby mode
5. Auto power off
6. Low power consumption
7. Low voltage indication

Contraindication

- Is not for continuous monitoring.
- Is not for those allergic to silicone, PVC, TPU, TPE, or ABS.
- Is not to be used if the skin tissue on/near the application site is damaged.
- Is not for use during cardiopulmonary resuscitation.
- Is not for use when the patient is hypovolemic.
- Is not for assessing the adequacy of ventilator support.
- Is not for detecting worsening lung function in patients on high concentrations of oxygen.

Intended Use


The Fingertip Pulse Oximeter with Wireless Connectivity is a handheld non-invasive device that is intended to be used for measuring oxygen saturation (SpO_2), and pulse rate through the user's finger. The product is suitable for spot-checking. Recommended for use by families at home, individuals in hospitals, oxygen bars, community healthcare clinics, or athletes (before or after physical exertion. Not recommended for use during workouts).

Product Accessories

1. Lanyard 2. User Manual 3. Carrying Case


To attach the lanyard use the following steps:

1. Insert the thinnest portion of the lanyard string through the small hole on the top side of the oximeter.
2. Hold open the section of the string that has been inserted through the hole to create an open loop.
3. Feed the large end of the lanyard through the smaller string loop.
4. Pull to tighten.

Operation Instructions

View of the front LCD


1. Install 2 x AAA batteries into the battery compartment correctly.
 a. Slide the battery door cover horizontally along the bottom of the pulse oximeter.
 b. Install two AAA batteries into the battery compartment. Match the polarities (plus (+) and minus (-) signs) in the compartment. If they are not matched, damage may be caused to the oximeter.
 c. Close the battery door.

Notes:
1. Please remove the batteries if the pulse oximeter will not be used for long periods of time.
2. Please replace the battery when the power indicator starts flickering.


2. Insert one finger into the rubber hole of the oximeter fully, making sure to place clamp over fingernail.
3. Press the button on the front panel to power the device ON.
 a. **Synchronous Time** will show while the device is turning on. This is attempting to synchronize with wireless connection. If you wish to use this device with the app, see section **MediLink App**. If you wish to proceed without connecting wait for this screen to pass.
4. Try to remain as still as possible while the oximeter is taking a measurement.
5. Read corresponding data from display screen.

When powered on, and in measurement mode, the device will automatically switch to one of four display modes depending on the direction the unit is held.
The displays are shown as follows:

If the oximeter does not detect a finger in the device, the screen will display "Finger Out".

MediLink App

For the most up to date information on what Android and iOS operating systems are compatible with our device applications please visit: biosmedical.com, the Apple Store, or the Play Store.

PLEASE NOTE: Android and iOS system and software updates may render the app obsolete on older devices.
This pulse oximeter and the **MediLink App** can work independently of each other. However, you will want to pair the oximeter and the App before use, as it will send the reading to the App once it is completed.

1. To begin using the free **MediLink App**, download and install onto your smartphone or tablet.

GET IT ON Google Play | Download on the App Store

2. Turn ON connection settings on your device.
3. Open the **MediLink App**. First time users, will need to register and create a password.
4. With the app open to the O_2 tab, select 'Add Pulse Oximeter'. (Figure 1)
5. Select your device from the list of Pulse Oximeters. (Figure 2).
6. Follow the instructions on your phone screen to pair the device.
 A pop up asking you 'Is this your pulse oximeter?' will appear. Select Yes. (Figure 3)
 ***Note:** Ensure that connectivity settings are enabled on your phone before attempting to pair.
7. Device has paired successfully when you can see the device listed at the bottom of your app home screen. (Figure 4)

Figure 1 **Figure 2** **Figure 3** **Figure 4**

NOTE: For more information on the App and how to use it, visit our website www.biosmedical.com or call our BIOS Medical Hotline: 1-866-536-2289

Maintenance and Storage

Turn off the device, and remove batteries before cleaning or disinfecting. Never immerse the device into liquid.

Cleaning/Disinfecting

Please use a soft cloth dampened with 75% alcohol to clean the silicone that touches the finger inside of the oximeter. Also clean the finger being tested with alcohol before and after each measurement.
DO NOT pour or spray liquids onto the oximeter, and do not allow any liquid to enter any openings in the device. Allow the oximeter to dry thoroughly before reuse.
CAUTION: Never use EtO or formaldehyde for disinfection.

Maintenance

- Check the main unit and all accessories periodically to ensure that there is no visible damage. It is recommended to inspect this device on a weekly basis.
- Replace the batteries in a timely manner when low indicator symbol is present.
- Clean the fingertip oximeter before and after each use.
- Remove the batteries from the device if the oximeter is not intended to be operated for an extended period.
- Keep in a dry place. Extreme moisture may affect oximeter lifetime and may cause damage.
- Dispose of batteries properly; follow any applicable local battery disposal laws.
- The fingertip pulse oximeter requires no routine calibration or maintenance other than replacement of batteries.

Specifications

Display Type	LCD Display	
SpO_2 *See notes 1-2	**Display Range:** 0% - 99% **Measurement Range:** 0% - 100% **Accuracy:** 70% - 100%: ±2%; 0% - 69%: unspecified **Resolution:** 1%	
Pulse Rate (PR) *See note 3	**Display Range:** 30bpm - 250bpm **Measurement Range:** 30bpm - 250bpm **Accuracy:** 30bpm - 99bpm, ±2bpm; 100bpm - 250bpm, ±2% **Resolution:** 1bpm	
Accuracy Under Low Perfusion *See note 4	**Low Perfusion:** 0.4% **SpO_2:** ±4% **PR:** 30bpm - 99bpm, ±2bpm; 100bpm - 250bpm, ±2%	
Light Interference	Under normal ambient light conditions, the SpO_2 deviation is ≤ 1%	
Pulse Intensity	Continuous bar graph display - The higher the display, the stronger the pulse	
Probe LED Specifications *See note 5	**Red Light:** Wavelength: ~660 nm; Optical Output Power: < 6.65 mW **Infrared Light:** Wavelength: ~905 nm; Optical Output Power: < 6.75 mW	
Wireless Connection Specifications	**Working Frequency:** 2402 MHz to 2480 MHz **Modulation Mode:** GFSK **Transmitting Power:** 0 dBm, +4 dBm **Receiving Sensitivity:** -93 dBm	
Operating Environment	**Temperature:** 10°C to 40°C (50°F to 104°F) **Relative Humidity:** ≤ 75% **Atmospheric Pressure:** 700 hPa to 1060 hPa	
Storage Environment	**Temperature:** -40°C to 60°C (-40°F to 140°F) **Relative Humidity:** ≤ 95% **Atmospheric Pressure:** 500 hPa to 1060 hPa	
Safety Class	Internally powered equipment, type BF applied part	
International Protection	IP22	
Working Voltage	DC 3.6 V	
Working Current	≤ 100 mA	
Power Supply	1.5V (AAA size) alkaline batteries x 2	
Operation Time	Can be continuously operated as long as 24 hours	
Dimensions	58(L) x 32(W) x 34(H) mm	
Weight	Approximately 52g (with batteries)	
Reference to Standards	IEC 60601-1:2012 IEC 60601-1-2:2014 EN 60601-1-6:2010 IEC 60601-1-11:2010 ISO 80601-2-61:2011 EN 13485:2012 EN ISO 14971:2012	EN ISO 14155:2011 EN 1041: 2008 EN 980:2008 IEC 62366:2007 IEC 62304:2006 EN ISO 10993-1:2009

NOTE 1: The claims of SpO_2 accuracy are supported by clinical study measurements taken over the full range. By artificially inducing a stable oxygen level to a range of 70% - 100% SaO_2, and then comparing against SpO_2 values collected by the secondary standard pulse oximeter equipment at the same time to form paired data, then used for the accuracy analysis.Twelve (12) healthy volunteers were used to collect data for the clinical report.

NOTE 2: Since pulse oximeter equipment measurements are statistically distributed, only about two-thirds (2/3) of the measurements can be expected to fall within ± Arms of the value measured by the co-oximeter.
NOTE 3: A patient simulator has been used to verify the pulse rate accuracy. It is stated as the root-mean-square difference between the PR measurements value and the value set by the simulator.
NOTE 4: A patient simulator has been used to verify the accuracy under conditions of low perfusion. SpO_2 and PR values are different due to low signal conditions. They are compared against the known SpO_2 and PR values of the input signal.
NOTE 5: Optical sensors such as the light-emitting components, will be affected by other medical devices with regard to the wavelength range. This information may be useful for clinicians carrying out optical treatment.

Possible Problems and Solutions

Problems	Possible reason	Solution
SpO_2 or PR values cannot be shown normally	• The finger is not properly inserted • The finger is shaking or the patient is moving • The device is not being used as outlined in the manual • Patient's SpO_2 value is too low to be measured	• Retry by inserting the finger • Try to calm down and stay still before measuring again • Read the manual and try measuring again • Try measuring again. If you are sure no problem exists with the product, please go to a hospital for timely diagnosis.
The device cannot be turned on	• The batteries are low, or completely drained • The batteries are installed incorrectly • The device is malfunctioning	• Replace with new batteries, and try measurement again • Ensure the battery polarities have been matched correctly • Contact Thermor Ltd. Customer Service
The display suddenly disappears	• The device has entered into energy saving mode • The batteries are low, or completely drained • The device is malfunctioning	• This is normal • Replace with new batteries, and try measurement again • Contact Thermor Ltd. Customer Service
Indication lamps are suddenly off	• The device has entered into energy saving mode • The battery power is too low to work	• This is normal • Replace the batteries

Symbol Definitions

Symbol	Definition	Symbol	Definition
	Follow instruction for use	Synchronous Time***	Synchronous time interface (Wireless equipment)
	Type BF applied part.	PR_{bpm}	Pulse Rate (bpm)
	Manufacturer's information	**%SpO_2**	Oxygen saturation
	Date of Manufacture		Low power indication
	Conformity to WEEE Directive		Temperature limitation
IP22	Protected against dripping water		Humidity limitation
Finger Out	Finger is not inserted in device, or not inserted correctly		Atmospheric pressure limitation
CE 0123	This device is compliant with Medical Device Directive 93/42/EEC	LOT	Lot No.

Notes:
- The illustrations used in this manual may differ slightly from the appearance of the actual product.
- The specifications are subject to change without prior notice.

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© 2025 Thermor Ltd. ALL RIGHTS RESERVED